

Haelvoet home collection

List with tested cleaning products:

- → Always follow the prescribed concentration
- → Always rinse with water and neutral soap after cleaning
- ightarrow Using cleaning agents can cause discoloration, there is no deterioration of the coating layer
- → Products marked with a " v " are listed on the VAH list of disinfectants

ТҮРЕ	Cleaner	Max. concentration
<u>Ammonia</u>		
	Sprint Multi E4c	100%
<u>Surfactants</u>		
	Reocid	5%
	1	'
ТҮРЕ	Surface disinfection	Max. concentration
<u>Alcohol</u>		
	Mikrozid [®] AF liquid v	6%
<u>Chlorine</u>	'	
<u></u>	Actichlor Plus	1000 ppm
	Chlor-Clean Tablets	1000 ppm
	Haz-Tabs	1000 ppm
	Sodiumhypochlorite	1000 ppm
Hydrogen peroxide		
	Hydrogen Peroxide	12%
	Incidin® OxyWipe S v	Wipes
<u>Oxygen</u>	1	'
	Perform ® v	6%
	Virkon Tablets	1%
Peracetic acid		
<u> </u>	Clinell Sporicidal wipes	Wipes
Quaternary		
<u> </u>	acryl-des® v	6%
	Incidin Foam v	Foam
	Surfa'Safe Premium v	Foam
	Umonium medical spray v	Spray
Quaternary & Alcohol	1	'
	Kohrsolin® extra v	6%